© Copyright MediRoyal Nordic AB, Sweden - 2021 - 25858DR - www.vajer.se

GB English – instructions

Mediroyal Pre-Shaped Hand Splints

General information - Read the instruction carefully. The product should be fitted by a medical professional like an occupational therapist, physical therapist, orthopedic technician, or a medical doctor. Patients with diabetes, circulatory problems or sensitive skin should consult their doctor before use. The splint is for single patient use and should not be reused on another patient than first fitted on.

Indications - The pre-shaped hand splints are designed to provide support and stability to the specific condition for each article. The splint might be an integrated part of the rehabilitation along with other splints or used as a stand alone treatment depending on the decision form the therapist.

TH1000 Reflex Inhibitory Position for patients with upper limb spasticity
TH1001 Cone splint for neutral position of the wrist and to inhibit spasticity over the fingers.
TH1002 Functional Positioning splint to support the wrist, and fingers.
TH1003 Neutral Positioning splint for neutral position of the wrist and abduction support to the thumb.
TH1004 Perforated Functional Positioning splint to support the wrist, and fingers.
TH1005 Resting Position splint to support the wrist and fingers.
TH1006 Functional Positioning splint to support the wrist, and fingers for children.

Contraindications - Patients that are sensitive to volar stimulation that can increase spasticity can be contraindicated for these splints. Patients with diabetes, material allergies, sensitive skin or poor circulation can also be contraindicated. The decision must be taken individually from case to case by the medical responsible staff.

Pre-application - Select the splint model depending on the diagnosis and the individual problem of the patient. To select the right size, use the measurements in the chart. Note that the product is available in left or right version.

Application and adjustment - The splint is made from a thermoplastic material that can be adjusted individually by using hot water (75–80 degrees C) or with a heat gun for spot heating. The splint is made from a low memory material that does not deform during heating. Never immerse the complete splint in hot water as the shape will be deformed. The splint might need further adjustments after a few days or use. Always check the function of the splint after a few days and observe pressure marks or edema. If these occur the splint must be adjusted.

Splint use and removal - The splint can be used both daytime and nighttime according to the decision from the medical responsible staff. To remove the splint, open all the straps.

Material content - The splint is made from Polycaprolactone CAS 24980-41-4

Caring instructions and precautions - We recommend that the splint should be washed once a week. Use mild soap and water. The splint can also be disinfected by using 70-85% rubbing alcohol. Patients with sensitive skin, diabetes or poor circulation should be observed. Patients with allergies should check the material content before use.

For more information visit www.mediroyal.se

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Preparation and adjustment - The splint is made from a thermoplastic material that can be adjusted individually by using hot water (75–80 degrees C) or with a heat gun for spot heating. The splint is made from a low memory material that does not deform during heating. Never immerse the complete splint in hot water as the shape will be deformed. The splint might need further adjustments after a few days or use. Always check the function of the splint after a few days and observe pressure marks or edema. If these occur the splint must be adjusted.

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Caring instructions and precautions - We recommend that the splint should be washed once a week. Use mild soap and water. The splint can also be disinfected by using 70-85% rubbing alcohol. Patients with sensitive skin, diabetes, or poor circulation should be observed. Patients with allergies should check the material content before use.

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