NL

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


Elbow brace Ulnaris · English

This brace has been designed to provide support and stability to the elbow. It can be used both to prevent flexion of the joint during the night as well as an active brace during the day to off-load the elbow region. The brace is fitted with both an aluminum stay and a thermoplastic piece to prevent flexion. You should only use one of them. The thermoplastic needs to be molded individually by a therapist or orthotist. The brace is also fitted with a pressure absorbing material over the elbow tip. The pocket is open in the end and the material can be taken out and adjusted individually.

- To select the right size use the measurement chart as a guide line. However if the brace should be used as a night splint you might want it to be a little bit bigger in size to not disturb circulation.
- Chose the flexion support that will work best; aluminum stay or thermoplastic.
- Aluminium –Bend the stay into the preferred position and slide it into the pocket. Make sure that it is completely inserted into the pocket.
- Thermoplastic –Use hot water, 75-80 degrees C to heat the plastic. It usually takes less than 1 minute to make it soft. If the patient has a sensitive skin, cover it with a sleeve before application. Apply the plastic onto the inside of the elbow and shape it into the preferred angle. Let it then cool and adjust the edges with a scissor. Insert the plastic into the pocket.
- Open the hook and loop straps on the brace.
- Slide the brace on to the elbow and close the straps. Don't tighten too hard.

If you have any allergic problems, check the material content on the backside of the label. Do not over- tighten the straps if you suffer from circulatory problems. If you are diabetic or suffer from skin problems, consult with your physician before using the brace. If your pain or problem gets worse, always consult with your therapist, physician or orthotist.

Material content: Lycra, Polyamide, Nylon, Spandex, PU-foam. This product is latex free.
Washing instruction: 40 degrees machine wash. See further instructions on product label.

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