

# MEDIROYAL OMOX FIX ARM SLING



**The Omox Arm Sling has been designed to provide an effective support and immobilization of the shoulder and arm after trauma.**

It's made of a stable 3D material that both absorbs pressure and provides a good ventilation. The neck strap is adjustable and has a pad to absorb pressure. The sling has also a thumb strap to provide support to the thumb and wrist. Around the waist is a foam strap that can be used to fixate the arm to the body to limit rotation of the shoulder joint.

#### INDICATIONS

Soft tissue trauma to the shoulder, upper- or lower arm when immobilization is needed. Can be used in combination with a cast to off-load the additional weight and to provide immobilization.

#### MR1102-9 MEDIROYAL OMOX FIX ARM SLING

The measurement is taken from the elbow tip to the MCP V joint.

Small	28-33 cm
Medium	33-40 cm
Large	40-51 cm