

# Genum Liga EZ

## Genum Liga EZ NeoTex Genum Liga EZ TriTex



EU4017 GENUM LIGA EZ NEOTEX  
EU4117 GENUM LIGA EZ TRITEX

Genum Liga EZ is a new knee orthosis with a full anterior opening. This makes the application easy since there is no need to pull the orthosis over the ankle. Fitted with a double hinged splint made of aluminum and steel that has been sewn in a NeoTex material for maximum protection. The splint is therefore completely encapsulated which gives a low profile and high function. The knee fold has a thin material made of NeoTex or TriTex to improve comfort. The popliteal area can be opened by using a pair of scissors. The material have to be cut evenly and no seams should be passed. The two wide Velcro straps have been fitted with additional closure straps on makes application easy. **INDICATIONS** Collateral Ligament Injuries, Meniscus Injuries, Osteoarthritis or other knee instabilities.

**X-SMALL** 29-32 centimeter  
**SMALL** 32-35 centimeter  
**MEDIUM** 35-38 centimeter

**LARGE** 38-42 centimeter  
**X-LARGE** 42-46 centimeter  
**XX-LARGE** 46-50 centimeter

 DeRoyal®

mediroyal 

Mediroyal Nordic AB  
Box 7052, 192 07 Sollentuna. Sweden  
Visiting address: Staffans väg 6B. Sollentuna  
Tel +46 8 506 766 00 Fax: +46 8 506 766 90  
[www.mediroyal.se](http://www.mediroyal.se) [info@mediroyal.se](mailto:info@mediroyal.se)