

QUEST AIR WALKER



Quest Air Walker

A walker with a full shell design, made of injection moulded plastic. Provide a cast-like fit to promote stability and patient comfort. The shell has an increased circumference to accommodate swelling. The soft-goods has an integrated air pump for patient and practitioner convenience. The wide base makes the walker stable but yet low profile. The sole has a rocker bottom to help promote a natural gait. The straps are adjustable and can be removed.

INDICATIONS Moderate to severe ankle sprains, Stable foot or ankle fractures.

1542 Quest Air Walker

| Size | Shoe size | Height | Circumference |
|---------|-----------|---------|---------------|
| Small | 35-39 | 37,5 cm | 38 cm |
| Medium | 39-42 | 40 cm | 48 cm |
| Large | 43-45 | 42,5 cm | 66 cm |
| X-Large | 45- | 42,5 cm | 68 cm |

