



mediroyal

Orthoses are our passion

CRX

CRX Epi Elbow

MR2810 · CRX Epi Elbow


The CRX Epi Elbow is made from a flat knitted compression textile with three medical graded compression zones. The highest compression zone is positioned over the joint, grading off into lower zone with less compression. The upper part of the support is made with a lower compression to not disturb range of motion.

On the inside of the brace are two silicone pads positioned, laterally and medially to provide compression over the extensor and flexor tendons in the lower arm. The silicone pads provide a combination of medical compression and a light to moderate support.

Indications: Lateral or Medial Epicondylitis.

Measurement chart

MR2810 CRX Epi Elbow

Size	XS	S	M	L	XL	XXL	XXXL	
cm	19–21	21–23	23–25	25–27	27–29	29–31	31–33	

NRX251 · NRX® Epi Elbow Brace 

© Mediroyal Nordic AB 2021

Mediroyal Nordic AB · Box 7052 · 192 07 Sollentuna · Sweden · Tel +46 8 506 766 00
info@mediroyal.se · facebook.com/Mediroyalnordic · www.mediroyal.se



NRX®

NRX® Epi Elbow Brace

NRX251 · NRX® Epi Elbow Brace

This is a new elbow brace for Epicondylitis made from NRX® material and Neptune™ with a modular strap and pad. The anatomical design makes it extremely functional for most activities. The lower part of the arm is made with a loop surface that allows the pad to be positioned exactly where is needed over the extensor muscles. The pad is made from 5 mm CR-foam with hook surfaces on both sides. The elastic strap is positioned on the pad and can then be adjusted individually for the right compression. The brace can also be used without the strap and pad once the epicondylitis problem has improved to insulate the body heat over the elbow joint.

The arm fold part of the brace has a Neptune™ material which is hydrodynamic, breathable, and quick drying. The upper and lower part of the brace has hook & loop adjustments.

Indications: Lateral or Medial epicondylitis.

Measurement chart

NRX251 · NRX® Epi Elbow Brace

Size	XS	S	M	L	XL	XXL	
cm	20–23	23–25	25–28	28–31	31–33	33–36	

mediroyal

Vajer | write | 2021 | -25906DR | ENG

MR2810 · CRX Epi Elbow