



mediroyal



Mediroyal Materials

Ventus™ Fleece Edging Tape and Padding

Our popular Ventus™ Fleece is now available as a self-adhesive edging tape and padding material. The Ventus™ fleece is made from spun polyester which provides a soft and comfortable surface that allows moisture transportation. The self-adhesive surface is allergy tested and has a water base which adheres to most textiles and thermoplastic materials.

Ventus™ Edging Tape 20 mm

EU451020 Ventus™ Edging Tape
20 mm width · 5 meter per roll
Black

Ventus™ Padding Material 30 cm

EU451300 Ventus™ Padding Material
30 cm width · 1 meter per roll
Black

Neptune™ Edging Tape and Padding



Mediroyal Materials

Neptune™ Edging Tape and Padding

The Mediroyal Neptune™ material is used in our popular Proxi® thumbs and is a hydrodynamic and functional material with high comfort. Now it is also available as self-adhesive edging tape and padding material. The material is made from nylon, polyester and spandex and has a high elasticity. It is also only 0,4 mm and extremely thin, making it possible to use in custom orthotics with very limited space. The self-adhesive surface is allergy tested and has a water based glue which adheres to most textiles and thermoplastic materials.

Neptune™ Edging Tape 20 mm

EU461020 Neptune™ Edging Tape
20 mm width · 5 meter per roll
Grey

Neptune™ Padding Material 30 cm

EU461300 Neptune™ Padding Material
30 cm width · 1 meter per roll
Grey



◀ **Ventus™ Fleece Edging Tape and Padding**