



DeRoyal
Sports Medicine



ANATOMICAL DESIGN
 HEAT AND COMPRESSION
STABILITY FLEXIBILITY
SUPPORTING

LIVING PROOF

DeRoyal by
mediroyal

no matter what...

No matter what kind of sports you perform there is always a risk for an injury. In the DeRoyal assortment you can find supports that can prevent an injury. We also have supports if you're already injured. The right type of support can reduce pain and inflammation so you can return to activity much faster.

SHOULDER



Shoulder support 8010 An anatomically designed shoulder support manufactured in NeoTex to provide heat, compression and stability to the injured shoulder. The support can be used for both left or right due to the unique strapping system. • Size XS-L



Shoulder support 8011 An anatomically designed shoulder support with additional straps to stabilize the shoulder. Manufactured in NeoTex to provide heat, compression and support. The support can be used for both left or right due to the unique strapping system. • Size XS-XL

ELBOW

Elbow support 8020 Anatomically designed to provide optimal fit and compression. The support is made of 3 mm NeoTex to provide heat, compression and support. Velcro closure for easy application. • Size XS-XXL



Elbow support Hyper-Ex 2625 Elbow support designed to prevent hyper-extension of the elbow. Made in 3 mm NeoTex together with stable hinges and X-strap to effectively prevent hyper-extension. SoftTex in the elbow crease for a high comfort. • Size XS-XXL



Elbow support Handball 8101 A support specially designed for handball. Made of 3 mm NeoTex and a shock absorbing pad covered with Kevlar® which provides extreme durability. • Size XS - XXL



WRIST






Wrist support with splint 8030 Made in one-piece design for the best fit and function. Velcro closure on the dorsal side along with a wide strap provides a good support. Supplied with a pre-shaped plastic splint. • Size S - XL Right and Left



Wrist strap Universal 8031 A flexible wrist strap that provides light support. Suitable for activities like tennis, golf and other sports that requires a free palm of the hand. The additional strap provides an adjustable stability. • Size Universal

BACK

Back support 8040 A light back support made of NeoTex to provide heat, compression and support. Velcro closure in the front for easy application. • Size S - XXL



Back support Stable 8041 A stable back support made of NeoTex and stabilizing splints in the back. An individual cross-strap provides compression and support. Velcro in the front for easy application. • Size S - XXL

...we're supporting you

THIGH



Thigh support 8060 Anatomical thigh support made of NeoTex to provide heat and compression. Velcro adjustment in the front. • Size S-XXL

KNEE

Knee Brace Basic

8050 Anatomically designed knee support for an optimal support and fit. Made of 5 mm NeoTex to provide heat and compression. The brace supports the knee and increases the circulation. • Size XS – XXL



Knee Brace with Open Patella

8051 Knee support with an open patella. Made of 5 mm NeoTex to provide heat and compression. The brace supports the knee and increases the circulation. • Size XS – XXL



Knee Brace with flexible stays


8052 Knee support made in 5 mm NeoTex to provide heat and compression. Flexible spiral stays on the side and a U-shaped patella pad provides additional support. Semi-elastic straps with velcro closure can be individually adjusted. • Size XS – XXL



Knee brace with stabilizing hinges

8053 Knee support designed to provide support to instable knees. Made of 5 mm NeoTex. The double-hinges are made of dural aluminium. They are placed in open pockets, enable them to be customized. Semi-elastic straps with velcro closure can be individually adjusted. • Size XS – XXL

Condro Band 8054 A small yet very effective brace. Helps to unload the patella tendon. Made of NeoTex Shark skin which grips the skin to prevent migration. • Size S – XL




Knee brace Handball 8100 Specially designed for handball. Made of 5 mm NeoTex and SoftTex in the knee fold. The pad is made of PE-foam covered with Kevlar® which provides extreme durability. • Size XS – XXL



Handball Knee brace Youngsters MR8102

Similar design as 8100 but adapted to active young players and made of 3 mm NeoTex. The pad is specially designed, made of PE foam, covered with Kevlar® for extreme durability. • Size 1 – 4

Volleyball Knee brace MR8004 Knee brace specially designed for Volleyball. Completely open in the knee fold to make it easy to flex the knee. Made of 5 mm NeoTex. The pad is made of PE-foam covered with Kevlar® which provides extreme durability. Sold in pairs. • Size S – XL



CALF


Calf support 8070 A anatomically designed calf support that provides heat and compression. Made of 3 mm NeoTex. Excellent for lower leg problems like shin-splints. Velcro closure in the front. • Size XS – XXL



SHORTS


Thermal pants

8110 Navy blue/8111 Black A true anatomical design that allows no restrictions in movement. Made of 1 mm NeoTex that provides both heat and compression. Excellent for groin, thigh or muscle tear problems. • Size XS – XL




Goal-keeper pants

8112 Same design as 8110 but with protective and shock absorbing pads for the hips. • Size S – XL



Compression pants

8115 The ultimate compression pants for light stabilization of the groin and thigh. Made of a thin and flexible Spandex Lycra • Size S – XL



ANKLE

Ankle support 8071

A soft and flexible ankle support made of 3 mm NeoTex. Provides light support and compression over the ankle and the achilles. Velcro closure in the front. • Size XS – XL



Ankle support with stabilizing strap 8072

Same model as 8071 but with an extra stabilizing strap to provide additional compression. Velcro closure. • Size XS – XL



Ankle brace 8075

Ankle brace designed for high impact activities. Made of ballistic nylon with low profile to provide optimal stability. Plastic stays for additional medial and lateral stability. The non-elastic stabilizing straps can be individually adjusted to lock the heel effectively to provide more stability. • Size XS – XL, Left/Right



Functional Ankle Brace 8073


A very stable ankle brace made in injection moulded plastic. Provides excellent stability and function both medially and laterally without disturbing dorsal or plantar flexion. Two elastic velcro straps provides an even compression over the ankle. The pads are made of PE foam. • Size XS – XL



SOLE

MediRoyal Impact

A highly functional assortment of shock absorbing insoles. Several different designs available like heel-spur wedges, shock absorbing insoles, arch supports and posting wedges. For more information, consult the MediRoyal Impact brochure.



TAPE

Kinesio tape Kinesio tape that provides correction to muscle tension, increases circulation and increases proprioception. To be able to use the full effect of the tape it's necessary to know how it works. For more information consult the MediRoyal Kinesio tape leaflet.



no matter what
...we're supporting you

Shoulder support

8010 Size	XS	S	M	L	–	–
8011 Size	XS	S	M	L	XL	–
Size in cm	24–28	28–32	32–35	35–38	38–41	–

Elbow supports

8020 Size	XS	S	M	L	XL	XXL
2625 Size	XS	S	M	L	XL	XXL
8101 Size	XS	S	M	L	XL	XXL
Size in cm	20–23	23–25	25–28	28–31	31–33	33–36

Wrist braces

8030 Size	–	S	M	L	XL	–
Size in cm	–	13–16	16–18	18–22	22–26	–
8031 Size	Universal					

Back supports

8040 Size	–	S	M	L	XL	XXL
8041 Size	–	S	M	L	XL	XXL
Size in cm	–	70–80	80–90	90–100	100–110	110–120

Thigh supports

8060 Size	–	S	M	L	XL	XXL
Size in cm	–	48–50	50–52	52–54	54–56	56–60

Knee braces

8050 Size	XS	S	M	L	XL	XXL
8051 Size	XS	S	M	L	XL	XXL
8052 Size	XS	S	M	L	XL	XXL
8053 Size	XS	S	M	L	XL	XXL
8054 Size	–	S	M	L	XL	–
8100 Size	XS	S	M	L	XL	XXL
8004 Size	–	S	M	L	XL	XXL
Size in cm	29–32	32–35	35–38	38–42	42–46	46–50
8102 Size	1	2	3	4	–	–
Size in cm	23–25	25–27	27–30	30–33	–	–

Calf support

8070 Size	XS	S	M	L	XL	XXL
Size in cm	34–35	35–36	36–37	37–38	38–39	39–40

Thermal and Compression pants / Shorts

8110 Size	XS	S	M	L	XL	–
8111 Size	XS	S	M	L	XL	–
8112 Size	–	S	M	L	XL	–
8115 Size	–	S	M	L	XL	–
Size in cm	60	70	80	90	100	–

Ankle supports

8071 Size	XS	S	M	L	XL	–
8072 Size	XS	S	M	L	XL	–
Size in cm	18–20	20–22	22–24	24–26	26–28	–
8073 Size	–	S	M	L	–	–
8075 Size	XS	S	M	L	XL	–
Shoe size	35–37	37–39	39–40	40–42	42–44	–

DeRoyal by

mediroyal

Mediroyal Nordic AB
Box 7052, Staffans väg 6B, 192 07 Sollentuna
Sweden

Tel: +46 8 506 766 00. Fax: +46 8 506 766 90

Website: www.mediroyal.se

E-mail: info@mediroyal.se