

TRANSITION TELE SCOPING POST-OP KNEE BRACE



A Post-Op knee brace with length-adjustable arms.

The length can be adjusted from 46-68 cm. Uses clip buckles for easy application and removal of the straps. The hinge can be set in 10 degree increments in both extension and flexion. Has also a drop-lock function that enables to lock the hinge in a fully extended **INDICATION** Post-operative stabilization and rehabilitation of knee injuries.

KB7000 Telescopic Post-Op

Size	Universal
Length	46-68 cm



mediroyal 

Mediroyal Nordic AB
Box 7052, 192 07 Sollentuna. Sweden
Visiting address: Staffans väg 6B. Sollentuna
Tel +46 8 506 766 00 Fax: +46 8 506 766 90

www.mediroyal.se
info@mediroyal.se
facebook.com/Mediroyalnordic